

Portavo 904 pH

Portable and robust analyzer for pH and ORP measurement

Up to 5,000 values can be recorded using the integrated data logger. Using the USB port and the Paraly SW 112 software, the logger data can be easily transferred to a PC for analysis.

Tailor-Made pH Calibration

Cal SOP

pH sensors can be tested using the SOP calibration procedure with up to 3 calibration points. The third buffer is used as a verification buffer. You can pick and choose the buffer set for each calibration point and also set the sequence. You can use your own buffer solutions or choose from a list of commercially available buffer sets, e.g., CaliMat, NIST and DIN. For the verification buffer, a maximum permitted deviation (delta pH) is entered.

Security Package Included

User management

The professional user management regulates access to the device and the sensor.

- Increased security for configuration, calibration and measurement data
- No unauthorized interventions during the operating cycle
- Up to 4 user profiles can be entered
- Different access rights can be established

Depending on user experience, the role profile can be selectively defined for configuring the device and sensor as well as for calibrating the sensor. The risk of changing settings inadvertently is clearly minimized in this way.

More Safety During Operation

Memosens sensors can be assigned directly to the Portavo 908 Multi.

As such, data saved in the sensor can be consulted, including

Sensor type

TAG

Group

Unique sensor-to-device assignment reduces potential errors. This ensures that only the right sensors are used for the selected measuring point.

По вопросам продаж и поддержки обращайтесь:

Архангельск (8182)63-90-72
 Астана +7(7172)727-132
 Астрахань (8512)99-46-04
 Барнаул (3852)73-04-60
 Белгород (4722)40-23-64
 Брянск (4832)59-03-52
 Владивосток (423)249-28-31
 Волгоград (844)278-03-48
 Вологда (8172)26-41-59
 Воронеж (473)204-51-73
 Екатеринбург (343)384-55-89
 Иваново (4932)77-34-06
 Ижевск (3412)26-03-58
 Иркутск (395) 279-98-46
 Киргизия (996)312-96-26-47

Казань (843)206-01-48
 Калининград (4012)72-03-81
 Калуга (4842)92-23-67
 Кемерово (3842)65-04-62
 Киров (8332)68-02-04
 Краснодар (861)203-40-90
 Красноярск (391)204-63-61
 Курск (4712)77-13-04
 Липецк (4742)52-20-81
 Магнитогорск (3519)55-03-13
 Москва (495)268-04-70
 Мурманск (8152)59-64-93
 Набережные Челны (8552)20-53-41
 Нижний Новгород (831)429-08-12
 Казахстан (772)734-952-31

Новокузнецк (3843)20-46-81
 Новосибирск (383)227-86-73
 Омск (3812)21-46-40
 Орел (4862)44-53-42
 Оренбург (3532)37-68-04
 Пенза (8412)22-31-16
 Пермь (342)205-81-47
 Ростов-на-Дону (863)308-18-15
 Рязань (4912)46-61-64
 Самара (846)206-03-16
 Санкт-Петербург (812)309-46-40
 Саратов (845)249-38-78
 Севастополь (8692)22-31-93
 Симферополь (3652)67-13-56
 Таджикистан (992)427-82-92-69

Смоленск (4812)29-41-54
 Сочи (862)225-72-31
 Ставрополь (8652)20-65-13
 Сургут (3462)77-98-35
 Тверь (4822)63-31-35
 Томск (3822)98-41-53
 Тула (4872)74-02-29
 Тюмень (3452)66-21-18
 Ульяновск (8422)24-23-59
 Уфа (347)229-48-12
 Хабаровск (4212)92-98-04
 Челябинск (351)202-03-61
 Череповец (8202)49-02-64
 Ярославль (4852)69-52-93



Facts and Features

- Memosens sensors or analog sensors for pH or ORP measurement can be used on one device.
- A sensor quiver protects the sensor from damage and drying out
- Robust housing with IP66/67, also for outdoor use
- Li-ion battery – charged directly via USB
- Data logger with 5,000 values
- Micro USB port and Paraly SW 112 software
- Tailor-made pH calibration Cal SOP
- User management for access control
- Sensor check for clear sensor-to-device assignment via sensor type, TAG or "Group"
- Temperature probe adjustment in the Memosens sensor (offset correction)

3-year
warranty!

MEMO SENS



Specifications

| | | |
|-------------------------------|---|--|
| pH/mV input (analog) | pH socket, DIN 19 262 (13/4 mm) | |
| | pH range | -2 ... 16 |
| | Decimal places*) | 2 or 3 |
| | Input resistance | $1 \times 10^{12} \Omega$ (0 ... 35 °C / 32 ... 86 °F) |
| | Input current | $1 \times 10^{-12} \text{ A}$ (at RT, doubles every 10 K) |
| | Measuring cycle | Approx. 1 s |
| | Measurement error ^{1,2,3)} | < 0.01 pH TC < 0.001 pH/K |
| | mV range | -1300 ... 1300 mV |
| | Measuring cycle | Approx. 1 s |
| | Measurement error ^{1,2,3)} | < 0.1% meas.val + 0.3 mV TC < 0.03 mV/K |
| Temperature input | 2 x 4 mm dia. for integrated or separate temperature detector | |
| | Measuring ranges | NTC 30 kΩ -20 ... 120 °C / -4 ... 248 °F Pt1000 -40 ... 250 °C / -40 ... 482 °F |
| | Measuring cycle | Approx. 1 s |
| | Measurement error ^{1,2,3)} | < 0.2 K (Tamb = 23 °C / 73.4 °F); TC < 25 ppm/K |
| Memosens pH, ISFET input | M8 socket, 4 pins, for Memosens lab cable | |
| | Display ranges ⁴⁾ | pH -2.000 ... 16.000 mV -2000 ... 2000 mV Temperature -50 ... 250 °C / -58 ... 482 °F |
| Sensor standardization*) | pH calibration | |
| Operating modes*) | Calmatic | Calibration with automatic buffer recognition |
| | Manual | Manual calibration with entry of individual buffer values |
| | Data entry | Data entry of zero and slope |
| | Cal-SOP (TAN option) | Software option SW-P001: Defining the pH buffers and the sequence of the calibration steps; defining the delta deviation for the verification buffer |
| | Temperature calibration (TAN option) | Software option SW-P002 for temperature probe adjustment in the Memosens sensor (offset correction) |
| Calmatic buffer sets*) | -01- Mettler-Toledo -02- Knick CaliMat -03- Ciba (94) -04- NIST technical -05- NIST standard -06- HACH -07- WTW techn. buffers -08- Hamilton -09- Reagecon -10- DIN 19267 -U1- (User) | 2.00/4.01/7.00/9.21 2.00/4.00/7.00/9.00/12.00 2.06/4.00/7.00/10.00 1.68/4.00/7.00/10.01/12.46 1.679/4.006/6.865/9.180 4.01/7.00/10.01/12.00 2.00/4.01/7.00/10.00 2.00/4.01/7.00/10.01/12.00 2.00/4.00/7.00/9.00/12.00 1.09/4.65/6.79/9.23/12.75 loadable via Paraly SW 112 |
| Permissible calibration range | Zero point | 6 ... 8 pH |
| | Slope | Approx. 74 ... 104 % |
| | Calibration timer*) | Interval 1 ... 99 days, can be switched off |
| | Sensoface | Provides information on the sensor condition |
| | Evaluation of | zero/slope, response, calibration interval |

Specifications

| | | | | | |
|---|---|---------------------------------|--|--|--|
| Memosens ORP input | M8 socket, 4 pins, for Memosens lab cable | | | | |
| Display ranges ⁴⁾ | mV | –2000 ... 2000 mV | | | |
| Sensor standardization* | Temperature | –50 ... 250 °C / –58 ... 482 °F | | | |
| ORP calibration (zero adjustment) | | | | | |
| Permissible cal. range | ΔmV (offset) | –700 ... 700 mV | | | |
| Temperature calibration (TAN option) | Software option SW-P002 for temperature probe adjustment in the Memosens sensor (offset correction) | | | | |
| Connections | | | | | |
| 2 x socket, 4 mm dia., for separate temperature probe | | | | | |
| 1 x M8 socket, 4 pins, for Memosens lab cable | | | | | |
| 1 x micro USB-B for data transmission to PC | | | | | |
| 1 x pH socket, to DIN 19262 | | | | | |
| Display | | | | | |
| LCD STN 7-segment display with 3 lines and icons | | | | | |
| Status indicators | For battery power level, logger | | | | |
| Notices | Hourglass | | | | |
| [on/off], [cal], [meas], [set], [▲], [▼], [STO], [RCL], [clock] | | | | | |
| Keypad | | | | | |
| Data logger | | | | | |
| 5,000 memory locations | | | | | |
| Recording | Manual, interval- or event-controlled | | | | |
| MemoLog calibration data logger (Memosens only) | | | | | |
| Up to 100 Memosens calibration records can be saved | | | | | |
| Recording | Directly retrievable via MemoSuite or Paraly SW 112 (USB) | | | | |
| Viewable on the display | Manufacturer, sensor type, serial no., zero, slope, calibration date | | | | |
| Communication | | | | | |
| USB 2.0 | | | | | |
| Profile | HID, driverless installation | | | | |
| Usage | Data exchange and configuration via Paraly SW 112 software | | | | |
| Diagnostic functions | | | | | |
| Sensor data (Memosens only) | Manufacturer, sensor type, serial number, operating time | | | | |
| Calibration data | Calibration date; zero and slope | | | | |
| Device self-test | Automatic memory test (FLASH, EEPROM, RAM) | | | | |
| Device data | Device type, software version, hardware version | | | | |
| Data retention | Parameters, calibration data > 10 years | | | | |
| EMC | | | | | |
| EN 61326-1 (General Requirements) | | | | | |
| Emitted interference | Class B (residential environment) | | | | |
| Immunity to interference | Industrial applications | | | | |
| EN 61326-2-3 | | | | | |
| RoHS conformity | According to directive 2011/65/EU | | | | |
| Power supply | | | | | |
| 4 x AA alkaline batteries | | | | | |
| 4 x AA NiMH battery or 1 x Li-ion battery, USB chargeable | | | | | |
| Operating time | Approx. 1000 h (alkaline) | | | | |
| Nominal operating conditions | | | | | |
| Ambient temperature | –10 ... 55 °C / 14 ... 131 °F | | | | |
| Transport/Storage temp. | –25 ... 70 °C / –13 ... 158 °F | | | | |
| Relative humidity | 0 ... 95 %, short-term condensing allowed | | | | |
| Housing | | | | | |
| Material | PA12 GF30 + TPE | | | | |
| Protection | IP 66/67 with pressure compensation | | | | |
| Dimensions | Approx. 132 x 156 x 30 mm / 5.2 x 6.14 x 1.18 inches | | | | |
| Weight | Approx. 500 g / 1.10 lbs | | | | |

* user-defined

1) at nominal operating conditions

2) ± 1 count

3) plus sensor error

4) ranges depending on Memosens sensor